## Test Accessories

Taber® Abraser & Shear / Scratch Tester Accessories

### Record test details on each specimen.

Specimen mounting cards (S-36 and S-36-1) provide a unique way to capture pertinent test data on each specimen. On the front of smooth card stock is a printed test report that enables critical information to be recorded. A pressure sensitive adhesive makes these cards ideal to mount flexible specimens. Plus, they prevent the formation of wrinkles caused by the rub-action of the abrading wheels. Mounting cards are offered in 4<sup>1</sup>/<sub>4</sub> inch (108mm) square or round formats.

## Mounting sheets provide non-slip attachment.

Coated on both sides with a pressure sensitive adhesive, specimen mounting sheets (S-37 and S-37-1) permit the testing of samples without a center hole. The adhesive sheet allows you to affix specimens directly to a specimen table that does not have the mounting screw. One side consists of a printed test report, permitting use of the sheet as a mounting card. Mounting sheets are offered in 4<sup>1</sup>/<sub>4</sub> inch (108mm) square or round formats.

### Specimen plates enable the testing of coatings.

Taber offers a selection of specimen plates to assist with the evaluation of coatings. Different substrates permit testing based on the usage of the product and include: wood, steel, aluminum, zinc, glass, and plastic.







# Test Accessories

Taber® Abraser & Shear / Scratch Tester Accessories

Specimen Mounting				
S-36	<b>Specimen Mounting Card, Square</b> $-4^{1/4^{u}}$ (108mm) square card coated on one side with a pressure sensitive adhesive, test record form printed on reverse side	50 pcs.		
S-36-1	Specimen Mounting Card, Round - same as S-36, but round	50 pcs.		
S-37	<b>Specimen Mounting Sheet, Square</b> $-4^{1/4"}$ (108mm) square card coated on both sides with a pressure sensitive adhesive, test record form printed on one side to permit using as mounting card	100 pcs.		
S-37-1	Specimen Mounting Sheet, Round – same as S-37, but round	100 pcs.		



Recommended storage conditions: temperature 23°C  $\pm 2^\circ C$  and relative humidity 50%  $\pm 5\%$ 

## **Specimen Plates**

S-16	<b>Specimen Plates, Steel</b> $- 4$ " (100mm) square, 20 gauge steel with a $1/4$ " (6.5mm) center hole	10 pcs.
S-17	<b>Specimen Plaques, Wood</b> – 4" (100mm) square, $\frac{1}{4}$ " (6.5mm) thick birch plywood with a $\frac{1}{4}$ " (6.5mm) center hole	10 pcs.
S-18	Specimen Plates, Aluminum – 4" (100mm) square, 19 gauge aluminum with a <sup>1</sup> /4" (6.5mm) center hole	10 pcs.
S-31	<b>NEMA Standardization Plates, Glass</b> – 4" (100mm) square, $\frac{1}{8}$ " (3mm) plate glass with a $\frac{1}{4}$ " (6.5mm) center hole	10 pcs.
S-34	<b>NEMA Standardization Plates, Zinc</b> $- 4$ " (100mm) square, $\frac{1}{32}$ " (0.8mm) zinc with a $\frac{1}{4}$ " (6.5mm) center hole	5 pcs.
S-38	<b>Specimen Plates, Plastic</b> – 4" (100mm) square, <sup>1</sup> / <sub>16</sub> " (1.5mm) Plexiglas with <sup>1</sup> / <sub>4</sub> " (6.5mm) center hole	10 pcs.





## Miscellaneous Rotary Abraser Test Accessories

S-11	<b>Refacing Disc</b> – 4" (100mm) diameter, 150 grit abrasive paper used to reface Taber Calibrase® abrading wheels	100 pcs.
ST-11	<b>Refacing Stone, Fine</b> $-4$ " (100mm) diameter x <sup>1</sup> /4" (6.5mm) thick, with <sup>1</sup> /4" (6.5mm) center hole, used to reface Taber Calibrase <sup>®</sup> CS-10F abrading wheels	1 each
140-109	<b>Counterweights</b> –for use with the Taber Rotary Abraser, reduces the arm load by 125g (175g counterweights also available)	2 each
	<b>Taber Clean-Up Hose</b> – attaches to the rear vacuum nozzle arm opening, facilitating instrument clean-up and cleaning of test specimens	



455 Bryant Street North Tonawanda, New York 14120 USA For more information: Call 1.800.333.5300 (U.S.) or 716.694.4000 (outside U.S.) Fax 716.694.1450 Email us at sales@taberindustries.com Visit us at www.taberindustries.com